

Assistive Hearing  
Solutions

**Wear&Hear**

PRODUCT PORTFOLIO



# OUR COMPANY



*When hearing fades, isolation grows. Hearing connects us to the rhythm of life, the voices of loved ones, the soaring beauty of music, nature and so much more. At some point we will all experience hearing degradation or know someone who is suffering from hearing loss - that's why we leveraged Alango's sound technologies to create the Wear & Hear product line.*

*Our hearing solutions package more than two decades of expertise in a way that is simple, easy to use and accessible for all. As engineers, we are privileged to apply our knowledge to connect people back to the heartbeat of experiences and literally, make life sound better.*

---

**Dr. Alexander Goldin,**  
**Founder & CEO**



# OUR PRODUCTS



# OUR PRODUCTS

## **BeHear NOW | ABH-101B** *Personalizable Bluetooth® hearing amplifier*

### When is it good?



Live  
Conversation



Mobile  
Calls



TV  
Listening



Music  
& Audio



Ambient  
Sounds


### Who is it for?

Active adults with mild-to-moderately severe hearing loss

### Main Features

- Customize to your hearing needs using the app
- Lightweight, with durable carrying case
- Four hearing presets: Indoor, Outdoor, Crowd and Live Music
- Rechargeable battery for up to 12 hours use

**BeHear®**  
Now



*Lightweight neck loop assistive hearing headset for active people*



**BeHear NOW | ABH-101B**  
*Personalizable Bluetooth® hearing amplifier*



# OUR PRODUCTS

**BeHear®**  
ACCESS

## **BeHear ACCESS | ABH-120** *Personalizable Bluetooth® hearing amplifier*

Bluetooth®



### When is it good?



Live  
Conversation



Mobile  
Calls



TV  
Listening



Music  
& Audio





Ambient  
Sounds

### Who is it for?

Older adults and those with dexterity issues  
who have mild-to-moderately severe hearing loss

### Main Features

- Customize to your hearing needs using the app
- Large buttons, charging cradle for ease of use
- Integrated telecoil receivers
- Six hearing presets: Indoor, Outdoor, Crowd, Live Music, T-coil Exclusive, T-coil Transparent
- Rechargeable battery for up to 13 hours use



*Neck loop assistive hearing headset  
with large controls and charging cradle  
(includes Telecoil receivers)*



**BeHear ACCESS | ABH-120**  
*Personalizable Bluetooth® hearing amplifier*

# OUR PRODUCTS

## **HearLink PLUS | BT5805A**

*Long-range transmitter that delivers optimized, personalized audio to BeHear headsets*

When is it good?



TV  
Listening



Music  
& Audio

Who is it for?

People who want to enhance their TV listening experience

### Main Features

- Low latency for optimal audio-video sync
- Long range connectivity (up to ~200 feet /~60 meters) for freedom of movement
- Seamless operation with BeHear headsets for remote play/pause control

**HearLink®**  
**PLUS**



*Long range, low latency,  
optical and analog wireless TV audio  
and music transmitter*





**HearLink PLUS | BT5805A**

*Long-range transmitter that delivers optimized, personalized audio to BeHear headsets*

# OUR PRODUCTS

## **BeHear PROXY | ABH-130**

*Personalizable Bluetooth® amplifier in a wireless neck speaker*

### When is it good?



Live  
Conversation



Mobile  
Calls



TV  
Listening



Music  
& Audio

### Who is it for?


People seeking a comfortable high-quality sound experience especially for TV watching

### Main Features

- Customize to your hearing needs using the app
- Ergonomic neck band rests on shoulders, leaving ears free
- Converts to personal amplifier and private listening device (using extractable earbuds)
- Suitable for use with or without hearing aids

# BeHear®

PROXY



*Comfortable wireless neck speaker  
with personalized amplification and  
retractable cables for listening privacy*



**BeHear PROXY | ABH-130**  
*Personalizable Bluetooth® amplifier in a wireless neck speaker*

# OUR PRODUCTS

## **BeHear BUSINESS | ABH-110** *Monaural personalizable amplifier*

### When is it good?



Live  
Conversation



Mobile  
Calls



Ambient  
Sounds

### Who is it for?

Professionals who desire a customized,  
discreet hearing solution

### Main Features

- Customize to your hearing needs using the app
- Hands-free phone calls
- Multi-function buttons for smartphone-independent operation
- User-defined directionality: omni, front, rear
- Operates in either ear



**BeHear®**  
**BUSINESS**



*Comfortable, lightweight, monaural headset  
with personalized amplification*






***BeHear BUSINESS | ABH-110***  
***Monaural personalizable amplifier***




# OUR APP

Customize your hearing experience for all BeHear products with the dedicated **W&H BeHear** app:


1 | *Personalize the headset using the app.*



2 | *Choose the hearing preset best matching the current environment.*



3 | *Fine tune what you hear in real time.*



## UNIQUE

BeHear headsets enable users to create a customized sound experience for live conversations, mobile calls, TV programs, ambient sound, etc.


# EXPRESS HEARING CHECK KIOSK

## Spark interest and engage new customers!

In less than 10 minutes customers can test their own hearing and compare the results with and without Wear & Hear personalized hearing amplification.

## Easy to install:

requires 12x10 inches of counter space,  
electricity, WiFi



## CONTACT US

---

Wear & Hear  
2 Etgar St., PO Box 62  
Tirat Carmel 3903213 ISRAEL  
Tel: +972 (0) 77-320-4102

[sales@WearAndHear.com](mailto:sales@WearAndHear.com)

**[www.WearAndHear.com](http://www.WearAndHear.com)**